


**BUET**



**EEE 426: Biomedical Signals,  
Instrumentation and Measurement  
Laboratory**



**An IoT Based Remote Health Monitoring  
System Using Integrated ECG  
Module and Pulse Oximeter**

# Group01: Lab Group G2



Farhan Hamid  
1706175



Sadia Afrose  
1706161



Shafim Bin Hassan  
1706172



Subah Karnine  
1706174



Joy Saha  
1706189



Md Al Amin  
Patwary  
1706127



Ishtiaque  
Ahmed  
1706190



# Aim Of This Project

- Creating a smart non-invasive heart monitoring system with integrated sensors
- Uploading the real time data to IoT cloud for remote access
- Finding heart rate from the acquired ECG signal and cross validating with a pulse oximeter



# Societal Impact

- Health care for people in living in remote areas without direct access to medical professionals
- A cheap home monitoring system for patients
- Easily portable alternative that will be convenient for use during different rescue operations in disaster zones







# Hardware Setup

Setup consisted:

- NodeMCU (ESP8266)
- AD8232
- 3-lead electrodes





# Final Hardware Setup






# Real Time Data Acquisition



# Identifying The Relevant Peaks from Acquired ECG signal





# Uploading BPM to ThingSpeak

Reading data from AD8232 for approx. 3 sec at 100sample/sec



Detrend and normalize the signal



Find the local maxima and identify the R peaks




Calculate the heart rate from mean R-R interval




Upload BPM to Thingspeak



Wait 30 seconds then repeat



# Calculating BPM



$$\text{Heart Rate} = \frac{Fs}{R-R \text{ interval}} * 60$$



# IoT Infrastructure



## ThingSpeak dashboard showing BPM graph and current value

### Channel Stats

Created: [22 days ago](#)

Last entry: [about 3 hours ago](#)

Entries: 9






# Sending ECG signal by mail



- Problem: Uploading the entire ECG signal to a IoT platform is not possible due to data loss and delay.
- Solution: Send the signal via mail with attached image of the signal and the related CSV file for further processing if desired.



# Sending Mail



# Mails sent under the two situations







# Uploading to Drive

Attached image file from the mail is automatically uploaded to Google Drive.

Pabbly connect is an intermediary mail parser that connects the two platforms



# Image uploaded to Google Drive



# Innovations and Findings

- Use of IoT cloud platforms in remote health monitoring
- Instant mailing for relevant information in case of abnormalities detected.
- Automatic information upload to Google drive





# Cost Sheet

Equipment	Cost(BDT)
AD3232 ECG MODULE	1800
Disposable ECG Electrodes pad	1200
ESP8266	345
Transport	500
Bread Board	70



# Future of the Project

- Provide greater reach to rural and underserved communities
- Enables frequent patient monitoring and data collection
- Don't have to drive to the doctor's office or clinic
- Help to prevent the spread of COVID-19, flu and other infectious diseases
- promote a healthier environment for everyone
- A family member can help by providing information, asking questions



# THANK YOU!!!

## Any Questions?

